

AMXDmax

Drainable Collection Bag




CB-10-2-RU-X



Features and Benefits

- Made from durable polyurethane
- Automatically collects urine as the Control Unit pumps it away from the Male Cup or Female Pad.
- Reusable for up to 20 flights when properly cleaned
- Liquid tight air filter and connectors
- 500mL over capacity safety factor
- Certified Safe to Fly



Catalog Number

CB-10-2-RU-12 - 1200ml Capacity
CB-10-2-RU-16 - 1600ml Capacity
CB-10-2-RU-20 - 2000ml Capacity

GSA

CCU-10 (RU-12)
CCU-10 (RU-16)
CCU-10 (RU-20)

Description

The Collection Chamber Unit is used with the Male Cup or Female Pad, Undergarment and system Control Unit to provide aircrew with in-flight bladder relief. It is drainable and reusable for up to 20 flights when properly cleaned.

Related Products

Removable Hose, AMXDmax, Through-Suit Connector, Collection Bag

Specifications

Capacity*	1200ml / 1600ml / 2000ml
Size	6.75"x0.5"x11.25"/14.25"/17.25"
Weight	1.9 oz / 2.6 oz / 3.3 oz
Operating Temp.	14° F to 122° F -10° C to 50° C
Storage Volume	1" x 1.5" x 6.5" 68.6mm x 86.4mm x 157.3mm
Storage Temp.	-4° F to 158° F -20° C to 70° C

* specify when ordering



Omni Medical Systems, Inc.

808 Hercules Drive, Colchester, VT 05446

☎: (802)497-2253 ☎: (802)497-3601

www.amxdmax.com