

# AMXDmax

## Control Starter Kit (Male)



### CSK-2 M



#### Contents

<b>Control Unit Kit</b>	SCU-10-2 M
<b>Collection Unit</b>	IMC-2-2* IMCR-2-2 w/ RHA-21B-2*
<b>Collection Bag</b>	CB-X-2* CB-10-2-RU-X* CB-10-2-LR-X*
<b>Undergarment</b>	MG-X*
<b>Disinfectant</b>	SDC-20

\* specify when ordering, see datasheet for options

#### Catalog Number

CSK-2 M

#### GSA

CSK (-2 Specify Male)

#### Description

The Control Starter Kit (Male) contains a complete AMXDmax system and enough supplies for proper training, providing fully automatic bladder relief. Certified Safe to Fly.

#### Contents

AMXDmax Male Control Unit Kit, a collection unit\*, a collection bag\*, a male undergarment\*, and cleaning fluid. Supply refills are available individually or through the Collection Supply Kit (CSK-B M)

#### Related Products

Collection Supply Kit

#### Specifications

<b>Weight</b>	3.5 oz
<b>Operating Temp.</b>	14° F to 122° F -10° C to 50° C
<b>Storage Volume</b>	4.5" x 3.5" x 3.75"
<b>Storage Temp.</b>	-4° F to 158° F -20° C to 70° C

\* specify when ordering



**Omni Medical Systems, Inc.**

808 Hercules Drive, Colchester, VT 05446

☎: (802)497-2253 ☎: (802)497-3601

www.amxdmax.com