

AMXDmax

Female Pad w/Built-in Hose



IFP-15-X-2



Features and Benefits

- Connects to the AMXDmax Female Control Unit for automatic operation
- Inflatable Ring Keeps Pad in Place
- Soft, Non-Woven Fabric Cover
- Urine sensor for automatic operation
- Worn discreetly under flight suit
- Non invasive, inflatable urethane pad
- Soft non-woven fabric cover for comfort
- Pad remains comfortably in place and ready to use for the duration of mission
- Flame retardant Nomex covered hoses
- Designed based on sanitary pad
- Certified Safe to Fly

Catalog Number

IFP-15-1-2-X*
IFP-15-2-2-X*

GSA

SCU-ABP-VA (15-1-X*)
SCU-ABP-VA (15-2-X*)

Description

The Pad w/Built-in Hose is an AMXDmax urine collection unit for females. It is made from flexible urethane and has a soft, non-woven fabric cover. The Pad inflates when connected to the AMXDmax Female control unit. The built in hose is covered with flame retardant Nomex fabric.

Related Products

AMXDmax, Through-Suit Connector, Collection Bag

Specifications

Hose Length	12" / 18" / 24"
Length	4.5" (144mm)
Width	3.4" (86mm)
Height	5.8" (157mm)
Weight	4.9 oz
Flow Rate	1.5LPM +/- .2LPM
Operating Temp.	14° F to 122° F -10° C to 50° C
Storage Volume	4.5" x 3.5" x 3.75" 144mm x 89mm x 95mm
Storage Temp.	-4° F to 158° F -20° C to 70° C

* replace "X" with hose size of 12, 18, or 24



Omni Medical Systems, Inc.

808 Hercules Drive, Colchester, VT 05446

☎: (802)497-2253 ☎: (802)497-3601

www.amxdmax.com