

AMXDmax

Female Pad w/Removable Hose



IFP-16-1-2

NSN: 1680-01-546-4426



Features and Benefits

- Detachable hose exits Top Center
- Uses Standard Removable hose or
- Uses TSC Removable Hose for use with ESC-35-2 Anti-Exposure Suit Adapter
- Connects to the AMXDmax Female Control Unit for automatic operation
- Soft, Non-Woven Fabric Cover
- Urine sensor for automatic operation
- Worn discreetly under flight suit
- Non invasive, inflatable urethane pad
- Pad remains comfortably in place and ready to use for the duration of mission
- Certified Safe to Fly

**Shown with RHA-21-B-2 F hose attached*

Catalog Number

IFP-16-1-2 (Size 1)**

***See Female Pad Undergarment Sizing Chart*

Description

The Inflatable Pad w/Removable hose is an AMXDmax urine collection unit (Pad) for females. It is made from flexible urethane and has a soft, non-woven fabric cover. The Pad is worn deflated during normal wear and inflates -creating a seal against the body when connected to the AMXDmax Female control unit to provide secure bladder relief. The removable hose exits vertically from the center of the pad.

Specifications

Length	4.5" (144mm)
Width	3.4" (86mm)
Height	5.8" (157mm)
Weight	3.1 oz
Flow Rate	1.5LPM +/- .2LPM
Operating Temp.	14° F to 122° F -10° C to 50° C
Storage Volume	4.5" x 3.5" x 3.75" 144mm x 89mm x 95mm
Storage Temp.	-4° F to 158° F -20° C to 70° C



Omni Medical Systems, Inc.

808 Hercules Drive, Colchester, VT 05446

☎: (802)497-2253 ☎: (802)497-3601

www.amxdmax.com