

# AMXDmax

In-Flight Bladder Relief  
Integrity Male Condom  
**IMC-2**



## HYDRATE AND DOMINATE THE SKY

**Eliminate Tactical  
Dehydration and  
maximize flight  
performance with  
the AMXDmax**

Tactical dehydration is dangerous.

Dehydration severely impacts G-tolerance, situational awareness, and decision making ability. It can affect anyone, even pilots with decades of flight experience. It only takes a second to lose your aircraft or your life.

The **AMXDmax** allows pilots to hydrate properly and relieve their bladder mid-flight, without interrupting the mission. Now fully-automatic and with a slimmer profile, the AMXDmax lets you reach your true potential.

**Omni Medical Systems, Inc.**  
808 Hercules Drive, Colchester, VT 05446  
☎: (802)497-2253 📠: (802)497-3601

[www.amxdmax.com](http://www.amxdmax.com)



# AMXDmax

Integrity Male Condom



## Disposable Urine Collection

### Designed for Comfort



The **IMC** (Integrity Male Condom) provides quick and easy bladder relief for pilots.

The condom sheath is latex-free and uses a hydro-colloid adhesive to keep it secure during flight. Sensors detect the presence of urine and automatically activate the **AMXDmax** control unit, evacuating the urine into a collection bag.

The **AMXDmax** is an effective solution for in-flight bladder relief and allows for longer missions without distraction.



The **Supply Kit** comes with:

*10x FLO Condoms*

*10x Disposable Collection Bags*

*1x Hose*

Use with removable hose.

Available in standard issue or TSC (Through Suit Connector) version

**Omni Medical Systems, Inc.**

808 Hercules Drive, Colchester, VT 05446

☎: (802)497-2253 📠: (802)497-3601

[www.amxdmax.com](http://www.amxdmax.com)

