

AMXDmax

Male Cup w/Built-in Hose



IMC-2-2-X



Features and Benefits

- Worn discreetly under flight suit
- Non invasive urethane cup
- Soft foam ring comfortably fits to body
- Enhanced Comfort & Flexibility
- Attaches to Collection Bag and AMXDmax Control Unit for fully automated use
- Reusable up to 20 flights with proper cleaning
- Fits like a flexible male athletic cup
- Cup stays comfortably in place and ready to Use for the duration of mission
- Hose lengths of 12", 18", and 24" available.
- Certified Safe to Fly

Catalog Number

IMC-2-2-12* - Hose Length 12"

IMC-2-2-18* - Hose Length 18"

IMC-2-2-24* - Hose Length 24"

GSA

SCU-RBP-VA (12" Hose)

SCU-RBP-VA (18" Hose)

SCU-RBP-VA (24" Hose)

Specifications

Length	2.7" / 68.6mm
Width	3.4" / 86.4mm
Height	5.8" / 157.3mm
Weight	4.64 oz
Flow Rate	1.5LPM +/- .2LPM
Operating Temp.	14° F to 122° F -10° C to 50° C
Storage Volume	2.7" x 3.4" x 5.6" 68.6mm x 86.4mm x 157.3mm
Storage Temp.	-4° F to 158° F -20° C to 70° C

* specify when ordering

Description

The Cup w/Built in Hose is an AMXDmax urine collection unit for males. It is made from medical-grade urethane and has a foam ring. The hose exit is at the top of the cup and connects to the AMXDmax Control Unit & Collection Chamber.

Related Products

Cup w/Removeable Hose



Omni Medical Systems, Inc.

808 Hercules Drive, Colchester, VT 05446

☎: (802)497-2253 ☎: (802)497-3601

www.amxdmax.com