

AMXDmax

Male Cup w/Removable Hose



IMCR-2-2



Features and Benefits

- Removable Hose
- Worn discreetly under flight suit
- Non invasive urethane cup
- Soft foam ring comfortably fits to body
- 20% slimmer profile than IMC-2-2
- Hose exits at the bottom of the cup
- 360° Swivel Hose Design
- Designed for Legacy & 5-Point Harness
- Compatible with the ESC-35-2
- Attaches to Collection Bag and AMXDmax Control Unit for fully automated use
- Reusable up to 20 flights with proper cleaning
- Cup stays comfortably in place and ready to Use for the duration of mission
- Certified Safe to Fly

Catalog Number

IMCR-2-2

GSA

SCU-RBP-B

Description

The Cup w/Removable Hose is an AMXDmax urine collection unit for males. It is made from medical-grade urethane and has a foam ring. The hose exit is at the base of the cup, and connects to a compatible hose (RHA-21B-2, RBU-40-2).

Related Products

Removable Hose, AMXDmax, Through-Suit Connector, Collection Bag

Specifications

Length	2.7" / 68.6mm
Width	3.4" / 86.4mm
Height	5.8" / 157.3mm
Weight	2.4 oz
Flow Rate	1.5LPM +/- .2LPM
Operating Temp.	14° F to 122° F -10° C to 50° C
Storage Volume	2.7" x 3.4" x 5.6" 68.6mm x 86.4mm x 157.3mm
Storage Temp.	-4° F to 158° F -20° C to 70° C



Omni Medical Systems, Inc.

808 Hercules Drive, Colchester, VT 05446

☎: (802)497-2253 ☎: (802)497-3601

www.amxdmax.com