


**Omni's Removable Hose
for Inside Connection to CCU-20 Suit Adapter Connector**

MFG Part: SCU-ABP-B

Product Name: Removable Hose RHA-21B-2-XX

AMXDmax Air Crew Ensemble Integration

- Single handed Connect and Disconnect
- Integrated Sensor interconnect
- Rotating Connection prevents Kinking
- 6", 12", 18" and 24" Hose lengths
- Connects to Male Cup / Female Pad
- Nomex Covered Hose



Single handed Connect and Disconnect

Omni's SCU-ABP-B Removable Hose is easily attached on the inside of the CCU-20 **Anti-Exposure/Chem Bio Suit Adapter** connecting to Omni's Male Cup or Female Pad.

- Specify (XX) 6", 12", 18" and 24" Hose length

The **Remote Collection Bag Umbilical Hose** is certified **Safe To Fly** on US Military Aircraft.

Available through GSA

Protected under US Patent **US 20080238093 A1** and other patents pending

- Weight : 18.4 g (0.65 oz)_
- Dimensions: 32 mm (1.25") Dia. X 10 mm (0.4") H
- Material: Flame Retardant Acetyl Copolymer
- Color: Connectors – Black, Nomex - Green

Omni Medical Systems, Inc, an Omni Measurement Systems, Inc. company

808 Hercules Drive, Colchester, VT 05446 802-497-2253 Fax: 802-497-3601

Sales@OmniMedicalSys.com