

Omni's Through Suit Adapter/Connector **Improved AMXDmax Aircrew Ensemble integration features**

MFG Part: CCU-20

Product Name: Anti-Exposure/Chem Bio Suit Adapter ESC-35-2

AMXDmax Air Crew Ensemble Integration

- Integrated sealed Check Valve
- Flight Suit
- CBRN Protection Ensemble
- Cold Water Immersion Garment
- Single handed Connect and Disconnect
- Quick Disconnect Emergency Release
- Integrated Sensor interconnect
- Retro-Fit any Air Crew Ensemble



Anti-Exposure/Chem Bio Suit Adapter Shown integrated with flight suit for F-35 Harness

The CCU-20 is for use with Omni CCU-21 Remote Collection Bag Umbilical Hose for connecting to Omni's Fully automatic, Hands Free AMXDmax™ In-Flight Bladder Relief System for ejection seat

Specify Air Crew Ensemble

The Anti-Exposure/Chem Bio Suit Adapter is certified **Safe To Fly** on US Military Aircraft.

Available through GSA

Protected under US Patent **US 20080238093 A1** and other patents pending

- Weight : 9.9 g (0.35 oz)
- Dimensions: 38 mm (1.5") Dia. x 12.7 mm (0.5") H
- Material: Flame Retardant Acetyl Copolymer
- Color: Black



Suit Adapter can be integrated easily into all Air Crew Ensembles. It is permanently Snap Locked into the fabric providing an Air Tight and/or Liquid Tight Environmental Seal.

Omni's SCU-ABP-B Removable Hose is easily attached on the inside of the Suit Adapter to connect to either the Male Cup SCU-RBP-B or the Female Pad SCU-ABP-VA

Omni Medical Systems, Inc, an Omni Measurement Systems, Inc. company
808 Hercules Drive, Colchester, VT 05446 802-497-2253 Fax: 802-497-3601
Sales@OmniMedicalSys.com