

AMXDmax

Disposable Pack (Holder)



SCU-ABP-2



Features and Benefits

- Battery Pack Holder for a 9v Lithium Battery
- Allows use of off-the-shelf battery alternative for the AMXDmax
- Certified Safe to Fly

Catalog Number

SCU-ABP-2

GSA

SCU-ABP (-2 AMXDmax)

Description

The Disposable Pack (Holder) allows the AMXDmax to use a disposable 9v Lithium battery.

Related Products

Rechargeable Battery Pack

Specifications

Length	2.38" (60mm)
Width	1.2" (30mm)
Height	0.71" (18mm)
Weight	0.9oz / 1.3oz
Operating Temp.	14° F to 122° F -10° C to 50° C
Storage Volume	2.38" x 1.2" x 0.71" 60mm x 30mm x 18mm
Storage Temp.	-4° F to 158° F -20° C to 70° C

* specify when ordering



Omni Medical Systems, Inc.

808 Hercules Drive, Colchester, VT 05446

☎: (802)497-2253 ☎: (802)497-3601

www.amxdmax.com